

PACSGEAR MEDIA WRITER

- ✓ Easy to configure, deploy and maintain
- ✓ Intuitive client software lets you work anywhere
- ✓ Burn or upload DICOM images, PDF and reports

Affordable, reliable DICOM CD/DVD burning.

Media Writer from Hyland Helthcare delivers affordable, reliable DICOM CD/DVD burning wherever images are created. From the desktop to the film room, Media Writer scales to meet all of your image distribution needs. Now, image-intensive departments such as radiology, cardiology, gastroenterology, dermatology, orthopedics and pathology have a distribution system that is easy to configure, deploy and maintain. An ideal solution for outlying departments, hospitals and imaging facilities, Media Writer is a low-maintenance, cost-effective and intuitive software/hardware solution that lets you work anywhere.

MEDIA WRITER FEATURES

- Directly burns or uploads DICOM CDs/DVDs from any PACS or modality user interface; intuitive client software lets you work anywhere.
- Uploads images to a variety of electronic image exchange services.
- Securely uploads images to a variety of cloud-based and point-to-point vendors.
- Safeguards your portable media with AES 256-bit encryption technology.
- For large studies, Media Writer spans multiple CDs or autodetects between CDs and DVDs.
- Includes Gear View Basic, an embedded DICOM viewer.
- Exports and anonymizes DICOM CDs/DVDs with one click.
- Creates customized splash screens and CD/DVD labels.
- Combines studies from multiple DICOM devices.
- Supports multipatient burning.
- Interface supports 12 languages.
- Supports IHE workflow.
- Media Writer includes additional features to scan documents, create electronic forms and import multimedia files — JPEGs, AVIs, MPEGs and DICOM CDs/DVDs.

The Hyland logo is positioned in the bottom left corner of the page. It consists of the word "Hyland" in a white, serif font, with a registered trademark symbol (®) to its upper right. The text is centered within a square background that features a vertical color gradient from green at the top to blue at the bottom.

**MEDIA WRITER
TECHNICAL SPECIFICATIONS**



Model	Media WriterD200	Media WriterD100
Burner/printer	Epson® Discproducer™	Epson® Discproducer™
Robotics	Yes, AcuGrip™ technology	Yes, AcuGrip™ technology
Typical usage	Hospitals, film rooms (up to 200 discs per day)	Departments, imaging centers, medical office buildings (up to 100 discs per day)
Burn/print time per disc	3-4 minutes, typical	3-4 minutes, typical
Disc capacity	Holds 100 discs (CD and/or DVD)	Holds 100 discs (CD and/or DVD)
CD/DVD drives	2 CD/DVD drives	1 CD/DVD drive
Recording speed	40x CD-R; 12x DVD±R; 8x DVD±R DL	40x CD-R; 12x DVD±R; 8x DVD±R DL
Media	Inkjet-printable 12cm	Inkjet-printable 12cm
Inkjet print method	MicroPiezo® technology	MicroPiezo® technology
Ink cartridge/type	Dye Ink/6 cartridge	Dye Ink/6 cartridge
Ink colors	Cyan, magenta, yellow, light cyan, light magenta, black	Cyan, magenta, yellow, light cyan, light magenta, black
Palette yield	16.7M colors, 1000 discs per six cartridges, typical	16.7M colors, 1000 discs per six cartridges, typical
Data interface	USB 2.0 high speed	USB 2.0 high speed
OS support	Windows XP (SP2), Windows 7 x32/x64	Windows XP (SP2), Windows 7 x32/x64
Electrical	100-120VAC/240VAC, 50/60Hz, 50W	100-120VAC/240VAC, 50/60Hz, 50W
Certifications	FCC Class B, CAN/CSA Class B, CE	FCC Class B, CAN/CSA Class B, CE
System requirements	1.4GHz processor, 10GB hard drive space, 1GB memory	1.4GHz processor, 10GB hard drive space, 1GB memory
Dimensions (HxWxD)	14" x 15" x 18" (348 x 377 x 465mm)	12.76" x 14.84" x 18.31" (324 x 377 x 465mm)
Weight	52.9 lbs. (24kg)	46.2 lbs. (21kg)


 EU Authorized Representative

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands


Manufacturer
Hyland LLC
 4309 Hacienda Drive, Suite 500
 Pleasanton, CA 94588 USA

Australian Sponsor:
 Emergo Australia
 201 Sussex Street, Level 20,
 Tower II Darling Park,
 Sydney, NSW 2000, Australia

For support, contact:
Tel: 1-844-535-1404
Web: <https://community.hyland.com>

Learn more at HylandHealthcare.com/EnterpriseImaging

Media Writer is not used for the purpose of diagnosis, cure, mitigation, prevention, or treatment of a disease or condition.

Part number: DOC-MW-BRO-EN-Rev A

©2020 Hyland Software, Inc. and its affiliates. All rights reserved. Hyland Healthcare is part of Hyland Software, Inc. and its affiliated companies. All Hyland product names are registered or unregistered trademarks of Hyland Software, Inc. or its affiliates in the United States and other countries.

Hyland®